

RHYNOSOFT

COMPUTARIZED RHINOMANOMETER

FAST AND ACCURATE



RHYNOSOFT

COMPUTARIZED RHINOMANOMETER

- Customized test
- Easy to install
- Patient database
- Automatic calculation of main parameters



TECHNICAL SPECIFICATIONS

CHANNELS	2 independent (pressure and flow)
SENSIBILITY	Pressure channel: 3Pa (0,3mm H ₂ O) Flow channel: 5cm ³ /S
APPLICATIONS	Automatic calculation of the air flow resistance. Automatic calculation of respiratory Automatic calculation of pressure and flow inhalation and exhalation peaks Automatic flow / pressure graphic generation
POWER SUPLAY	110V / 220V
PC CONNECTION	USB port
MINIMAL REQUIREMENTS	Pentium II 400 MHz, 64MB RAM, 1G of free disk space
OPERATING SYSTEM	Windows®98SE, Windows® XP, Windows® ME

ECLERIS RESERVES THE RIGHT TO MODIFY THIS PRODUCT WITHOUT PRIOR NOTICE

Rev.06.11.18