

# COLPOSCOPE HD

SERIES C-100

- Outstanding image quality • Perfect HD 1280X720p color definition • Optical zoom
- LED high power illumination • Great maneuverability



# COLPOSCOPE HD

## SERIES C-100

### Optional:

- Monitor stand pole
- Endodigi software capturing system.
- JPG output



### TECHNICAL SPECIFICATIONS

<b>OUTPUTS</b>	2x HDMI
<b>MAGNIFICATION</b>	Progressive magnification. Optical zoom up to 37X and digital zoom up to 70X on a 32" monitor.
<b>FOCUS</b>	Automatic and manual
<b>RESOLUTION</b>	1280 x 720P
<b>WORKING DISTANCE</b>	240-300 mm
<b>FILTER</b>	Optical green
<b>LIGHT SOURCE</b>	LED 6000K integrated on head
<b>LED life</b>	50000 hs.
<b>VOLTAGE</b>	100-240 volts.
<b>FREQUENCY</b>	50/60hz

- OSD (on screen display) with indication of real magnification.
- Zoom and focus operated by 4 commands joystick.
- Still function.
- WB (white balance) push bottom control.
- Perfect color discrimination due to CMOS HD technology.
- Floor stand system with 4 antistatic wheels with break.
- Ratched up down pole with double arm for easier and steady positioning.

ECLERIS RESERVES THE RIGHT TO MODIFY THIS PRODUCT WITHOUT PRIOR NOTICE