

LAP25

LAPAROSCOPIC INSUFFLATOR

- Heated CO2 output
- New user friendly display
- Full automatic self test



LAP25

LAPAROSCOPIC INSUFFLATOR

The LAP25 insufflator achieves a quick, stable and precise intra-abdominal pressure.

New user friendly display panel indicates the following parameters:

- Preset abdominal pressure.
- Real abdominal pressure.
- Real flow.
- Total gas consumption.

Extreme accuracy is insured through electronic and mechanical pressure control devices.

Visual and audio alarms immediately alert the user of any irregularity.

The LAP25 performs a self-test every time it is turned on.

The user can select both the pressure level and the flow rate before the surgery while the system is in standby mode.

The LAP25 has a built-in CO₂ heating device for patient safety and comfort.

3 automatic flow levels.

TECHNICAL SPECIFICATIONS

MAXIMUM GAS FLOW	25 l/min. continuous
PRESSURE RANGE	1-30 mmHg
CONSUMPTION	65 W
POWER SUPPLY	110/240 VAC
FREQUENCY	50-60 Hz
CO ₂	Heating System

ECLERIS RESERVES THE RIGHT TO MODIFY THIS PRODUCT WITHOUT PRIOR NOTICE