

PORTALED

THE TRUE PORTABLE LIGHT SOURCE



ECLERIS 
INNOVATION VALUE

PORTALED

THE TRUE PORTABLE LIGHT SOURCE

Specially designed for endoscopy.

✓ Available connectors for:

- Rigid Scopes Storz and Wolf type.
- Fiberoptic cables Storz and Olympus type.



✓ White powerful LED light 100% true colors (CRI: 0,9 +/-3%).

✓ Battery always fully charged by Ecleris pad charger.

✓ Stainless steel connections, durable.

✓ Belt holder.

✓ Carrying case.



Made in USA

ECLERIS 
INNOVATION VALUE

info@ecleris.com | www.ecleris.com

SCAN ME

