

# EVERLUX

SUPER BRIGHT **LED LIGHT SOURCE**

80/W 50,000 HOUR LED LIFE



BRIGHT, WHITE AND POWERFUL

# EVERLUX

SUPER BRIGHT LED LIGHT SOURCE

Interchangeable fiber optic connection.  
Compatible with main brands



## TECHNICAL SPECIFICATIONS

LED	80 W
Light color temperature	6000° k
Power consumption	175 watts
Class	I
Type	BF
LED life	50.000 hours approx.
Over temp. protection	Automatic LED system power-off
Dimmer	Electronic
Fiber optic connector	Universal, compatible with major brands
Dimensions	202 X 290 X 95mm (width x depth x height)
Weight	2.5 Kg.
Power supply	100 - 240 V
Frequency	50 / 60Hz