

ECLERIS 

PROCAM FULL HD


FULL HIGH DEFINITION ENDOSCOPY CAMERA

1920 X 1080



PROCAM FULL HD

- ERGONOMICALLY DESIGNED CAMERA HEAD (COOL AND COMPACT).
- AUTOMATIC WHITE BALANCE W/MEMORY FUNCTION.
- PRESET PARAMETERS FOR DIFFERENT MEDICAL APPLICATIONS.
- 3 FUNCTION KEYS AT THE CAMERA HEAD FOR CONTROLLING ZOOM BRIGHTNESS AND SNAPSHOT.
- EXCELLENT COLOR REPRODUCTION.
- SHARP IMAGES.
- SUPERIOR DEPTH PERCEPTION.
- USB PORT FOR RECORDING.



6 PRE-LOADED APPLICATIONS

Sinuscopy & Otoscopy, Laryngoscopy, Fiberscope, Laparoscopy, Cystoscopy & Hysteroscopy/ Arthroscopy



TECHNICAL SPECIFICATIONS

SERIES	PROCAM HD
	1920 X 1080
VOLTAGE	100-240 VCA 50/60 Hz
POWER	25W
SENSOR CHIP SIZE	1/3" Sony CMOS sensor chip
OUTPUT VIDEO	1920 x 1080p 60fps
OUTPUT PORTS	HDMI, DVI, VGA, SDI, CVBS
IMAGE DISPLAY RATIO	FL 16:9
IMAGE SHUTTER SPEED	1/60fl-fl1/15360 sec.
SNR	>= 50dB
MINIMUM ILLUMINATION	1LUXflF1.2
CONTROLS ON CAMERA HEAD	Zoom, Brightness, Picture
CONTROLS ON CCU	White balance, Speciality, Picture, Video
PICTURE FORMAT	1920 x 1080
VIDEO FORMAT	MP4

PROCAM FULL HD



info@ecleris.com | www.ecleris.com



Ecleris is a Trademark of Ecleris USA DBA of EUSA Global LLC.
11801 NW 100th. Rd. Suite 17. Medley FL 33178.
Phone +1 786 483 7490