

PROCAM HD

FULL HIGH DEFINITION ENDOSCOPY CAMERA

1920 X 1080P

- Ergonomically designed camera head (cool and compact).
- Compact and portable CCU
- Automatic white balance w/memory function.
- 2 function keys at the camera head for controlling zoom and brightness.
- Excellent color reproduction.
- Sharp images.
- Superior depth perception.



ECLERIS RESERVES THE RIGHT TO MODIFY THIS PRODUCT WITHOUT PRIOR NOTICE

Rev. 27.09.22



PROCAM HD

FULL HIGH DEFINITION ENDOSCOPY CAMERA

1920 X 1080P

Voltage	120-240 VCA 50/60 Hz
Power	20W
Sensor chip size	1/3" Sony Cmos Sensor Chip
Output Video	1920x1080p 60fps
Output Ports	2 x HDMI
Image display ratio	f1 16:9
Image Shutter_Speed	1/60fl ~fl ~1/15360 sec.
SNR	>= 50dB
Minimum_Illumination	1 LUXflF1.2
Controls in camera head	zoom. brightness
Controls on CCU	White balance



ECLERIS RESERVES THE RIGHT TO MODIFY THIS PRODUCT WITHOUT PRIOR NOTICE

Rev. 27.09.22

