

ECLERIS 

NEW 2K STROBOSCOPY

SYSTEM



HIGH POWER LED STROBE
WITH NEW GENERATION
VIDEOCAMERA TECHNOLOGY

www.ecleris.com

NEW 2K

STROBOSCOPY SYSTEM

EXCELLENT IMAGE QUALITY

The Ecleris LED strobe includes a 2064x 1544 p resolution camera, more than 3mp definition and 12 Bits x pixel providing extreme high depth colors and contrast, razor-sharp images of the vibrating vocal folds during endoscopic exams. The LED strobe light technology allows for the clear visualization of vocal fold edges.

ENHANCED ILLUMINATION

High Power LED and high performance optical system, provide great illumination for both stroboscopy and endoscopy images without the need of periodic bulb replacement.

FASTER, EASIER EXAMS

Endodigi Software together with Footswitch control for key functions allows physicians to perform faster and easier exams. Online remote support is available in order to help with service issues.



NEW 2K

STROBOSCOPY SYSTEM

COMPLETE SYSTEM

- Super powerful LED STROBE
- 2K Camera
- Recording
- Software



STROBOLED

Light	LED 90-W 6500 k
Light modes	Continuous + Automatic Sync with voice 60 to 1000 Hz Sync range
Strobe Modes	Locked and 2 speeds on slow motion
Microphone	Pneumatic + Air Ambient
Camera	3 MP CMOS (2064x1544 - 12 bits) with Global Shutter and Sync. Digital USB 3.0
Visual Interface	Endodigi Software on Windows PC

NEW 2K

STROBOSCOPY SYSTEM

- High power LED Stroboscope with both continuous and strobe light in a single unit.
- Automatic transition during fonation.
- Automatic detection of fundamental frequency shown on digital display.
- Automatic and manual frequency sync.
- Footswitch with dual function both automatic and locked modes.
- Universal connection for all kind of Fiber Optic cables.
- Light intensity control.
- Pneumatic and air microphone.
- Integrated Computer with 2K Videocamera and Patient management system.
- Touch screen front panel with voice waveform display.
- Video and Pictures recording and advanced video/picture edition.
- Complete Patient database.
- Slow/Fast motion function.
- Stroboscopy Cart with flexible and rigid scope holders.

LED 50.000
HOURS LIFE
NO BULB
REPLACEMENT



NEW GENERATION STROBOLED



info@ecleris.com | www.ecleris.com



Ecleris is a Trademark of Ecleris USA DBA of EUSA Global LLC.
11801 NW 100th. Rd. Suite 17. Medley FL 33178.
Phone +1 786 483 7490