

# EVERLUX UNIVERSAL

- Universal fiberoptic connection



- Portable and compact light source, with LED technology, 80 watts, 50,000 hours life LED.
- Cold and white, high intensity light, 6.000 K color temperature.
- Universal connection. For all type of FO cables.
- Electronic dimmer.

## TECHNICAL SPECIFICATIONS

LED	80 W
Light color temperature	6000° k
Power consumption	175 watts
Class	I
Type	BF
LED life	50.000 hours approx.
Over temp. protection	Automatic LED system power-off
Dimmer	Electronic
Fiber optic connector	Universal, compatible with major brands
Dimensions	202 X 290 X 95mm (width x depth x height)
Weight	2.5 Kg.
Power supply	100 - 240 V
Frequency	50 / 60Hz



[www.ecleris.com](http://www.ecleris.com)

**ECLERIS**   
INNOVATION VALUE